# GeneAce SYBR™ qPCR Mix α Low ROX

## I. Description

GeneAce SYBR™ qPCR Mix α Low ROX is a convenient 2X premix for real-time qPCR using SYBR™ Green I dye. This product was optimized for instruments that require a low ROX reference signal for fluorescence normalization (Compatible instruments: ABI 7500, Stratagene Mx3000P / 3005 / 4000, etc.).

This product provides high specificity and reproducibility and allows for reducing non-specific amplification by using optimized buffer components and Hot-Start Gene Taq NT, a chemically modified Taq DNA Polymerase.

## II. Storage

-20°C, protected from light.

This product should be protected from light. After thawing, the master mix may be stored at 4°C for up to one month or returned to -20°C for long term storage.

## III. Components

**300 rxns**

2 x GeneAce SYBR™ qPCR Mix α Low ROX *1)	1.5 ml x 5
---	------------

\*1) It contains Hot-Start Gene Taq NT, dNTP Mixture, Mg<sup>2+</sup>, SYBR™ Green I, ROX reference dye, stabilizer.

## IV. Notes

- **The PCR must start with an initial 10 minutes incubation at 95°C to activate the chemically modified hot-start Taq DNA polymerase.**
- **We recommend two-step cycling protocols. (See V. Typical PCR protocol)**
- **This product cannot be used for runs in “Fast”<sup>\*)2</sup> mode on ABI real time instruments.**

\*2) 95°C 20 sec → (95°C 3 sec , 60°C 30 sec) x cycle numbers

- Always ensure that the product has been fully thawed and mixed before use. Gently mix the mixtures without creating bubbles.
- Uracil-N-Glycosylase treatment is unusable on this product.

## V. Typical PCR Protocol

Choose an appropriate total reaction volume, depending on the instrument used. (e.g., 50 µl or 25 µl)

Component	Volume	Volume
2x GeneAce SYBR qPCR Mix α Low ROX	25.0 µl	12.5 µl
25 µM each Primers	1.0 µl	0.5 µl
Template	5.0 µl	2.5 µl
d.d.H <sub>2</sub> O	up to	50.0 µl 25.0 µl

Recommended cyclor conditions

**95°C 10 min. Enzyme activation <sup>\*)3</sup>**

**95°C 30 sec. ) 45 cycles**

**60°C 1 min. )**

Melting curve analysis

Optimum PCR condition depends on primer sequence and template DNA and so on.

\*3) Ensure that the cycling program includes the DNA polymerase activation step. Initial 10 minutes incubation at 95°C is required.

### For Licensee Products sold in the Research Field:

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

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