

---

---

# ISOSPIN Plasmid

Manual (Ver. 05)

---

---

Code No. 318-07991

NIPPON GENE CO., LTD.

## I Description \_\_\_\_\_

The ISOSPIN Plasmid is an easy and rapid method for the purification of highly pure plasmid DNA from 1-5 mL overnight culture of *E. coli*.

## II Kit components \_\_\_\_\_

Component	(100 preps)	Note
IS1 Buffer	30 ml x 1	(see next page)
IS2 Buffer	30 ml x 1	
IS3 Buffer	40 ml x 1	
ISPW Buffer	60 ml x 1	
ISW Buffer	100 ml x 1	(includes ethanol)*
ISE Buffer	10 ml x 1	10 mM Tris-HCl (pH 8.5)
RNase A (100 mg/ml)	60 µl x 1	(see next page)
Spin Column (a Spin Column and Collection Tube set)	50 sets x 2	

\* Keep the buffer bottle tightly closed after use.

## III Storage conditions \_\_\_\_\_

All the kit components can be stored at room temperature (15-25°C). RNase A is stable for longer periods when stored at 4°C or -20°C. After adding the RNase A to the IS1 Buffer, store the IS1 Buffer [+ RNase A] at 2-10°C.

## IV Precautions ---

- The kit is intended research use only.

## V Protocol ---

### **Material not supplied**

- Micropipette
- Pipette tips
- 1.5 mL microcentrifuge tubes
- Microcentrifuge

### **Before You Begin** ---

#### **Buffer Preparation**

- Add all the RNase A (100 mg/ml) to the IS1 Buffer, and mix well.

Store the IS1 Buffer [+ RNase A] at 2-10°C.

*Please mark the buffer bottle to indicate whether RNase A has been added or not.*

#### **Prepare the Starting Materials**

Pick a fresh single colony from an LB agar plate with appropriate antibiotic, and inoculate in liquid LB medium containing the same selective antibiotic. Incubate bacterial cultures at 37°C for 12-16 hours. Do not incubate the cultures for longer than 16 hours because the cells will begin to lyse.

You can now isolate your plasmid DNA from the bacterial culture:


- Prepare 1-5 ml *E.coli* culture for High-copy plasmids, such as pUC.
- Prepare 1-10 ml *E.coli* culture for Low-copy plasmids, such as pBR.

#### **Checklist**


- Please check whether the RNase A had been added to the IS1 Buffer bottle.
  - If precipitates appear in the IS2 Buffer, please incubate at 37°C, inverting periodically to dissolve.
-

## Plasmid DNA Extraction Protocol


1. Prepare **1-5 ml overnight (12-16 hours) culture of *E.coli***, as described in the previous page. Add the culture to a 1.5 mL microcentrifuge tube. Centrifuge at 10,000 x g for 5 min at room temperature. Discard the supernatant and repeat Step 1 to collect more cells.



2. Add 250  $\mu$ l of **IS1 Buffer [+ RNase A]** to the pellet of bacterial cells, and mix by vortex or pipet. Resuspend completely the pellet until there are no visible clumps.
3. Add 250  $\mu$ l of **IS2 Buffer** and gently invert the tube 4-6 times to mix until the solution is clear.  
Note) Do not allow this step for more than 5 min. Do not Vortex to avoid shearing of genomic DNA.
4. Add 350  $\mu$ l of **IS3 Buffer** and gently invert the tube 4-6 times to mix well.  
Note) Do not vortex or vigorous shake.
5. Centrifuge at 12,000 x g, for 10 min at room temperature.




6. Prepare a **Spin Column** (a Spin Column and Collection Tube set). Apply all the **supernatant** to the Spin Column.  
Note) Carefully transfer ~800  $\mu$ l supernatant to Spin Column without disturbing a white pellet form.
7. Centrifuge at 12,000 x g, for 1min at room temperature.
8. Discard the flow-through. Place the Spin Column back into the same Collection Tube.




9. Add 500  $\mu$ l of **ISPW Buffer** to the Spin Column.
10. Centrifuge at 12,000 x g, for 1min at room temperature.
11. Discard the flow-through. Place the Spin Column back into the same Collection Tube.






12. Apply 750  $\mu$ l of **ISW Buffer** to the Spin Column.
13. Centrifuge at 12,000 x g, for 1min at room temperature.
14. Discard the flow-through. Place the Spin Column back into the same Collection Tube.



15. Centrifuge at 12,000 x g, for 1min at room temperature to remove residual liquid.



16. Transfer the Spin Column to a new 1.5 mL microcentrifuge tube.
17. Add 50  $\mu$ l of **ISE Buffer** to the center of the Spin Column. Incubate 3 min at room temperature.  
Note) Nuclease-free water or TE (pH8.0) can also be used to elute the DNA in place of ISE Buffer (10 mM Tris-HCl, pH 8.5). We recommend using TE (pH8.0) containing EDTA if you store the eluted DNA for long term storage.
18. Centrifuge at 12,000 x g, for 1min at room temperature.
19. Recover your purified plasmid DNA in the microcentrifuge tube.  
Note) The purified plasmid DNA can be used directly or stored at  $-20^{\circ}\text{C}$  for long term storage.

## Simplified Protocol

1-5 ml overnight (12-16 hours) culture of *E.coli*

↻ Centrifuge at 10,000 x g, for 5 min at room temperature.  
Discard the supernatant.

Bacterial cells pellet

← Add 250  $\mu$ l of IS1 Buffer [+ RNase A] to the pellet, and mix by vortex or pipet.  
Resuspend completely the pellet until there are no visible clumps.

← Add 250  $\mu$ l of IS2 Buffer, and invert the tube 4-6 times to mix until the solution is clear.  
*Do not vortex. Do not leave for more than 5 min.*

← Add 350  $\mu$ l of IS3 Buffer and gently invert the tube 4-6 times to mix well. *Do not vortex.*

↻ Centrifuge at 12,000 x g, for 10 min at room temperature.

Apply all the supernatant to a Spin Column. Be careful to pipette the supernatant only.

↻ Centrifuge at 12,000 x g, for 1min at room temperature.  
Discard the flow-through.

← Add 500  $\mu$ l of ISPW Buffer to the Spin Column.

↻ Centrifuge at 12,000 x g, for 1min at room temperature.  
Discard the flow-through.

← Apply 750  $\mu$ l of ISW Buffer to the Spin Column.

↻ Centrifuge at 12,000 x g, for 1min at room temperature.  
Discard the flow-through.

↻ Centrifuge at 12,000 x g, for 1min at room temperature to remove residual liquid.

Transfer the Spin Column to a new 1.5 mL microcentrifuge tube.

← Add 50  $\mu$ l of ISE Buffer to the center of the Spin Column.  
Incubate 3 min at room temperature.

↻ Centrifuge at 12,000 x g, for 1min at room temperature.

Recover your purified plasmid DNA in the microcentrifuge tube.

## VI Troubleshooting

---

### **Low DNA Yield**

Ensure the ISE Buffer is applied to the center of the Spin Column so that elution is efficient. For Low-copy plasmid, use 1-10 ml of an overnight *E.coli* culture as starting material.

### **Contaminated with genomic DNA**

Do not vortex or leave after adding the IS2 Buffer. Please read instructions carefully.

### **Low DNA Performance**

If ethanol has been carried-over, spin for 5 min, instead of 1 min in Step 15.

## VII Data

---

### **Technical Information**

Binding Capacity:	up to 20 µg
Column Volume:	900 µl
Elution Volume:	50 µl
Plasmid Size:	up to 20 kbp

The information in the descriptions of the products may be changed without prior notification.

**NIPPON GENE CO., LTD.**

If you have any questions, please contact us by web form.

<http://www.nippongene.com/>