EasyPrep DNA Extraction Reagent

Code No. 391-4310

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Code No. 391-4310

I Product Description

This product is DNA extraction reagent. DNA can be extracted by adding this product to the animal tissues and by heat treating and by centrifuging. It takes only 15 minutes. The extracted DNA can be used as a template DNA for Loop-mediated Isothermal Amplification (LAMP).

[Precaution]

1. This product cannot be used for medical practice and clinical diagnosis.
2. This product should be handled in accordance with the descriptions in the manual. We are not responsible for problems caused if this product is not handled in accordance with the manual.
3. This product should be handled only by persons having a basic knowledge about reagents.
4. Expire date is printed on the label.
5. Follow the instructions from your region to discard this product.
6. Read SDS before use. SDS can be downloaded from our web site (<http://nippongene-analysis.com>).
7. LAMP is patented by Eiken Chemical Company.

II Contents

For 15 mL (50 reactions)

	Contents	Volume	Qty
I	EasyPrep DNA Extraction Reagent	15 mL	1
II	Manual	—	1

Storage

- ◆ This product can be stored at room temperature but should be stored at 4°C for long term storage. Do not freeze because a precipitate might be generated by freezing-melting.

III Necessities in addition to this product

- Micropipette and tip
- Centrifuge
- Micro tube and rack
- Heat block
- Pestle
- Vortex
- Crush ice
- Sterilized water

IV Protocol

1. Extraction of DNA from tissue

- 1) Weigh 50-300 mg[†] of muscle to 1.5 mL tube.
- 2) Add 300 µL of EasyPrep DNA Extraction Reagent and grind and mix with pestle, and then vortex well.
- 3) Incubate for 10 min at 95°C.
- 4) Incubate for 1 min on ice.
- 5) Centrifuge at 15,000×g for 5 min at room temperature.
- 6) Transfer supernatant into a new 1.5 mL tube.

Option: If extracted DNA inhibits nucleic amplification reaction, the supernatant should be diluted to 2-10 times by sterilized water.

[†] The volume of EasyPrep DNA Extraction Reagent can be changed depending on the kinds of the tissues.