



# UniversAll™ Extraction Buffer II

Cat. No. FYU004-5ML, SYU004-001

Storage: -20 °C

The UniversAll™ Extraction Buffer II is designed to enable DNA amplification directly from a wide range of biological samples without DNA purification. Simply incubate the tissue in the extraction buffer for 10 minutes at 98 °C and then use 1 µl of lysate for PCR. No Proteinase K treatment is necessary.

## DNA extraction

**Important Note: White precipitate is crucial for the extraction. Please make sure the white precipitate is added when add the buffer to the samples.**

- (1) Please Mix the buffer well before adding into samples. Add 50 µl of the UniversAll™ Extraction Buffer II to each tissue sample (**-1 mm<sup>3</sup> or 1 mg or 2 µl**) in a microcentrifuge tube. For samples like 0.1X serum, saliva and sputum, up to 50 µl of the samples can be extracted using 50 µl of the buffer.
- (2) Vortex 5 sec and centrifuge briefly. For the solid sample, make sure the sample block is submerged in the buffer.
- (3) Heat at 95-98°C for 10 min.
- (4) Vortex 5 sec and centrifuge briefly.
- (5) Use 1-2 µl of DNA extract for PCR amplification or qPCR analysis.
- (6) For those difficult samples that contain paraffin, phenolic compounds, heavy metals or some unknown inhibitory metabolites, a 10~100X serial dilution of the lysate is recommended before PCR amplification. The dilution can be done simply using PCR-grade water.

