GM Soybean (LLS) Detection Plasmid Set -ColE1/TE-

Code No. 315-07381 (1 Set)

Storage:

Store at -20°C, or at -80°C for long term storage.

Component:

The kit contains five concentrations of the standard plasmid and ColE1 DNA/TE as non-template control.

ColE1 DNA/TE Solution	(ColE1 DNA 5ng/μL)	90 μL×2
LLSoy-ColE1/TE-20	(20 copies/2.5μL)	90 μL×2
LLSoy-ColE1/TE-125	(125 copies/2.5μL)	90 μL×2
LLSoy-ColE1/TE-1.5k	(1,500 copies/2.5μL)	90 μL×2
LLSoy-ColE1/TE-20k	(20,000 copies/2.5μL)	90 μL×2
LLSoy-ColE1/TE-250k	(250,000 copies/2.5μL)	90 μL×2

Usage example:

Use 2.5 μL of the plasmid in each 25 μL reaction (e.g. Add 2.5 μL of this plasmid to 22.5 μL of PCR mixture.)

Note:

One package of the kit is controlled the quality as one product. We cannot guarantee quantitativity of the standard curve if reagents from other packages are used together.

Storage precautions:

- Low concentrations of DNA can be unstable if stored for long periods. We recommend storing the kit at -80°C for long-term storage.
- Avoid repeated freeze-thaw cycles more than 6 times.
- Do not aliquot the plasmid for storage, as there is a risk of contamination.
- All reagents should be stored under the same conditions; do not separate it away from the product package.

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
2. Takabatake R, Akiyama H, Sakata K, Onishi M, Koiba T, Futo S, Minegishi Y, Teshima R, Mano J, Furui S, Kitta K.: *Shokuhin Eiseigaku Zasshi*, **52** (2), 100 (2011)


Figure: Structure of GM Soybean (LLS) Detection Positive Control Plasmid

This product is a linearized plasmid diluted with ColE1 DNA/TE Solution.

Primer pair and Probe name:

GM Soybean(LLS) Detection LLS Oligonucleotide Set (Code No. 312-07391)

GM Soybean(RRS) Detection Le1(for endogenous gene) Oligonucleotide Set (Code No. 319-05601)

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