GeneAce Probe qPCR Mix II

I. Description

GeneAce Probe qPCR Mix II is a convenient 2X premix for real-time qPCR using sequence-specific fluorogenic probes.

This product enables reliable results across a wide dynamic range and allows for reducing non-specific amplification by using optimized buffer components and Hot-Start Gene *Taq* NT, a chemically modified *Taq* DNA Polymerase. This product contains dUTP and therefore carryover contamination is prevented by adding UNG (This product does not contain UNG).

This product contains a unique passive reference dye that is compatible across a variety of instrument platforms.

II. Storage

-20°C, protected from light.

This product should be protected from light. After thawing, the master mix may be stored at 4°C for up to one month or returned to -20°C for long term storage.

III. Components

	50 rxns	200 rxns
2 x GeneAce Probe qPCR Mix II ^{*1)}	1.25 ml x 1	1.25 ml x 4

*1) It contains Hot-Start Gene *Taq* NT, dNTP Mixture with dUTP, Mg²⁺, Passive reference dye, stabilizer.

IV. Notes

- **The PCR must start with an initial 10 minutes incubation at 95°C to activate the chemically modified hot-start *Taq* DNA polymerase.**
- **We recommend two-step cycling protocols. (See V. Typical PCR protocol)**
- **This product cannot be used for runs in “Fast”^{*2)} mode on ABI real time instruments.**

*2) 95°C 20 sec→(95°C 3 sec, 60°C 30 sec) x cycle numbers

- Always ensure that the product has been fully thawed and mixed before use. Gently mix the mixtures without creating bubbles.

V. Typical PCR Protocol

Choose an appropriate total reaction volume, depending on the instrument used. (e.g., 50 µl or 25 µl)

Component	Volume	Volume
2 x GeneAce Probe qPCR Mix II	25.0 µl	12.5 µl
25 µM each Primers	1.0 µl	0.5 µl
10 µM TaqMan [®] Probe ^{*3)}	1.0 µl	0.5 µl
Template	5.0 µl	2.5 µl
d.d.H ₂ O	up to	50.0 µl
		25.0 µl

Recommended cycler conditions ^{*4), *5)}

95°C 10 min. Enzyme activation ^{*6)}
95°C 30 sec.
60°C 1 min.) 45 cycles

Optimum PCR condition depends on primer sequence and template DNA and so on.

*3) TaqMan[®] is registered trademarks of Roche Molecular Systems, Inc.

*4) Depending on the primer sequence, the annealing time can be decreased to 15 sec, and the extension time can be decreased to 30 sec (optional).

*5) For carry over prevention, add UNG in reaction mix and incubate prior to qPCR reactions (optional).

*6) Ensure that the cycling program includes the DNA polymerase activation step. Initial 10 minutes incubation at 95°C is required.

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