# GeneAce SYBR™ qPCR Mix II

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# GeneAce SYBR™ qPCR Mix II

## I. Description

GeneAce SYBR™ qPCR II is a convenient 2X premix for real-time qPCR using SYBR™ Green I dye.

This product provides high specificity and reproducibility and allows for reducing non-specific amplification by using optimized buffer components and Hot-Start Gene Taq NT, a chemically modified Taq DNA Polymerase.

This product contains dUTP and therefore carryover contamination is prevented by adding UNG (This product does not contain UNG).

This product contains a unique passive reference dye that is compatible across a variety of instrument platforms.

## II. Storage

-20°C, protected from light.

This product should be protected from light. After thawing, the master mix may be stored at 4°C for up to one month or returned to -20°C for long term storage.

## III. Components

	125 rxns	500 rxns
2 x GeneAce SYBR™ qPCR Mix II <sup>(*)</sup>	1.25 ml x 1	1.25 ml x 4

<sup>(\*)</sup>1) It contains Hot-Start Gene Taq NT, dATP, dCTP, dGTP, dUTP, Mg<sup>2+</sup>, SYBR™ Green I, Passive Reference Dye, Stabilizers.

## IV. Notes

- **The PCR must start with an initial 10 minutes incubation at 95°C to activate the chemically modified hot-start Taq DNA polymerase.**
- **We recommend two-step cycling protocols. (See V. Typical PCR protocol)**
- Always ensure that the product has been fully thawed and mixed before use. Gently mix the mixtures without creating bubbles.


## V. Typical PCR Protocol

Choose an appropriate total reaction volume, depending on the instrument used.

Component	Volume [20 µl rxn]
2x GeneAce SYBR™ qPCR Mix II <sup>(*)</sup>	10.0 µl
25 µM each Primers	0.4 µl
Template	2.0 µl
d.d.H <sub>2</sub> O	up to 20.0 µl

<sup>(\*)</sup>2) UNG is not included in the GeneAce SYBR™ qPCR Mix II. UNG can be purchased separately from Nippon Gene (Code No. 317-09041). If optional treatment with UNG is required, add UNG to the reaction mixture (see the UNG manual).

Recommended thermal cycling conditions. <sup>(\*)</sup>3), <sup>(\*)</sup>4)

<sup>(\*)</sup>3) Optimal PCR condition depends on primer sequence and template DNA and so on.

<sup>(\*)</sup>4) Recommended ramp rates are 2.63-4.4°C/sec (up ramp) and 2.2-3.5°C/sec (down ramp). When using Thermo Fisher Scientific thermal cycler, we recommend selecting the "Fast mode" ramp speed instead of Standard mode. Ensure that the cycling program includes the initial 10 min 95°C step required to activate the DNA polymerase. Also, the extension time should be 1 min at 60°C.

<sup>(\*)</sup>5) For PCR fragments <200 bp, this denature time is 5 sec. For PCR fragments 200-500 bp, the denature time is 15 sec. For optimal results, PCR products should be between 80-180 bp and should not exceed 500 bp.

### For Licensee Products sold in the Research Field:

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.