GeneAce Turbo RT Kit with DNase I

Code No. 315-09841

保存: Store at -20°C

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Storage: Store at -20°C

Product Description

GeneAce Turbo RT Kit with DNase I enables rapid, high-yield cDNA synthesis for real-time PCR after genomic DNA removal from RNA samples.

Features

- Removes residual genomic DNA by DNase I treatment.
- Reverse transcription completed in 10 minutes at 50°C.

Contents

For 100 reactions (20 µL reaction volume)

Component	Volume
① DNase I (RNase free)	20 µL × 1
② 14× DN Turbo Buffer	100 µL × 1
③ 5× RT Turbo Buffer	400 µL × 1
④ RT Turbo Primer Mix	100 µL × 1
⑤ Turbo Enzyme Mix	100 µL × 1
⑥ ddWater (RNase free)	1.4 mL × 1

*RT Turbo Primer Mix contains Random Hexamers and Oligo(dT) Primer premixed at an optimized ratio. Depending on the application, Oligo(dT) Primer alone, Random Primer alone, or a gene-specific primer may be used instead.

*Turbo Enzyme Mix contains a modified M-MLV Reverse Transcriptase and RNase Inhibitor.

Precautions

- Mix reagents gently by inversion. Avoid vortexing or excessive agitation that may cause foaming.
- Keep all reagents on ice during reaction setup. Immediately return reagents to storage at -20°C after use.
- Real-time PCR reagents are not included in this kit and must be supplied separately.
- Use this kit before the expiration date.

Protocol

Step 1. Genomic DNA Removal

Prepare the following reaction mixture on ice.

① DNase I (RNase free)	0.2 µL
② 14× DN Turbo Buffer	1.0 µL
RNA Template (up to 1 µg)	up to 12.8 µL
⑥ ddWater (RNase free)	to 14.0 µL

Mix gently and incubate at 37°C for 5 minutes. After incubation, place the reaction on ice.

Note: If DNase I Treatment Is Omitted

14× DN Turbo Buffer and RNA template must still be added to the Step 2 reaction mixture.

Step 2. Reverse Transcription

Prepare the following reaction mixture on ice.

③ 5× RT Turbo Buffer	4.0 µL
④ RT Turbo Primer Mix	1.0 µL
⑤ Turbo Enzyme Mix	1.0 µL
Total	6.0 µL

Because 5× RT Turbo Buffer stops DNase I activity, promptly add the Step 1 reaction mixture to the mixture prepared above after the genomic DNA removal reaction.

Mixture prepared above	6.0 µL
Step 1 Reaction Mixture	14.0 µL
Total	20.0 µL

Mix and incubate under the following conditions:

25°C	5 min	Primer annealing
50°C	10 min	Reverse transcription
85°C	5 min	Enzyme inactivation
4°C	Hold	Storage

The synthesized cDNA may be used immediately for real-time PCR or stored at -20°C.

Note: The recommended cDNA input is ≤10% of the total PCR reaction volume. Optimization may be required depending on the PCR system used.

For inquiries, please use the contact form on our website.

<https://www.nippongene.com/>