


# Lambda DNA



酵素名	切断数	切断箇所
<i>Acc</i> I	9	0-100%
<i>Ava</i> I	8	0-100%
<i>Apa</i> I	1	0-100%
<i>Bal</i> I	18	0-100%
<i>Bam</i> H I	5	0-100%
<i>Bgl</i> II	6	0-100%
<i>Bss</i> H II	6	0-100%
<i>Bst</i> E II	13	0-100%
<i>Dra</i> I	13	0-100%
<i>Eco</i> O109 I	3	0-100%
<i>Eco</i> R I	5	0-100%
<i>Eco</i> T38 I	7	0-100%
<i>Hpa</i> I	14	0-100%
<i>Hind</i> III	7	0-100%
<i>Kpn</i> I	2	0-100%
<i>Mlu</i> I	7	0-100%
<i>Nar</i> I	1	0-100%
<i>Nco</i> I	4	0-100%
<i>Nde</i> I	7	0-100%
<i>Nhe</i> I	1	0-100%
<i>Nru</i> I	5	0-100%
<i>Nsi</i> I	14	0-100%
<i>Pvu</i> II	15	0-100%
<i>Sac</i> I	2	0-100%
<i>Sac</i> II	4	0-100%
<i>Sal</i> I	2	0-100%
<i>Sca</i> I	5	0-100%
<i>Sma</i> I	3	0-100%
<i>Sph</i>	6	0-100%
<i>Ssp</i> I	20	0-100%
<i>Stu</i> I	6	0-100%
<i>Sty</i> I	10	0-100%
<i>Xba</i> I	1	0-100%
<i>Xho</i> I	1	0-100%