



GM Soybean (RRS) Detection Le1n02 (for endogenous gene) Oligonucleotide

Code No. 313-05501 (60 µL)

Code No. 319-05503 (180 µL)

Storage:

Store at -20°C

Component:

Primer pair (25 µM each) 1 tube

Detection of:

Le1 (endogenous gene)

Primer name:

Le1-n02 [Le1n 02-5' & Le1n 02-3']

Amplicon length:

118 bp

Concentration:

This product contains both forward and reverse primers in a tube. The final concentration for each primer is 25 µM.

Reference:

1. "JAS Analytical Handbook, Manual of Inspection and Analysis for GM-Foods"; in Japanese, by IAI Center for Food Quality, Labeling, and Consumer Services of Japan.
2. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.
3. Kuribara, H., Shindo, Y., Matsuoka, T., Takubo, K., Futo, S., Aoki, N., Hirao, T., Akiyama, H., Goda, Y., Toyoda, M. and Hino, A.: *J. AOAC Int.*, **85** (5), 1077 (2002)
4. Kodama, T., Kuribara, H., Minegishi, Y., Futo, S., Watai, M., Sawada, C., Watanabe, T., Akiyama, H., Maitani, T., Teshima, R., Furui, S., Hino, A. and Kitta, K. : *J. AOAC Int.*, **92** (1), 223 (2009)

ダイズ内在性 DNA Le1n02 オリゴヌクレオチド

Code No. 313-05501 (60 µL)

Code No. 319-05503 (180 µL)

保存:

Store at -20°C

内容:

Primer pair (25 µM each) 1本

検出:

Le1 (ダイズ内在性遺伝子)

プローブ&プライマー名:

Le1-n02 [Le1n 02-5' & Le1n 02-3']

増幅長:

118 bp

濃度:

25 µM each

本品は、フォワードプライマーとリバースプライマーの混合溶液です。

参考文献:

1. 農林水産消費技術センター「JAS 分析試験ハンドブック 遺伝子組換え食品検査・分析マニュアル」
2. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」
3. Kuribara, H., Shindo, Y., Matsuoka, T., Takubo, K., Futo, S., Aoki, N., Hirao, T., Akiyama, H., Goda, Y., Toyoda, M. and Hino, A.: *J. AOAC Int.*, **85** (5), 1077 (2002)
4. Kodama, T., Kuribara, H., Minegishi, Y., Futo, S., Watai, M., Sawada, C., Watanabe, T., Akiyama, H., Maitani, T., Teshima, R., Furui, S., Hino, A. and Kitta, K.: *J. AOAC Int.*, **92** (1), 223 (2009)

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。