



GM Papaya Detection 55-1(Event-specific) Oligonucleotide Set

Code No. 310-07711 (1 Set)

Storage:

Store at -20°C and protect from light.

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 60 µL

Detection of:

PRSV-cp (55-1)

Probe & primer name:

PRSV-cp F & PRSV-cp R

PRSV-cp P (FAM/TAMRA)

Reference:

1. "Detection Methods of approved GM Foods"; in Japanese, by Consumer Affairs Agency.

GM パパイヤ系統別 DNA 55-1(Event-specific) オリゴヌクレオチドセット

Code No. 310-07711 (1 Set)

保存:

Store at -20°C(遮光)

内容:

Primer pair (25 μM each) 180 μL

Probe (10 μM) 60 μL

検出:

PRSV-cp (55-1)

プローブ&プライマー名:

PRSV-cp F & PRSV-cp R

PRSV-cp P (FAM/TAMRA)

参考文献:

1. 消費者庁「安全性審査済みの組換え DNA 技術応用食品の検査方法」

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。