



# GM Potato Detection J3 Oligonucleotide Set

---

**Code No. 318-08591** (1 Set)

**Storage:**

Store at -20°C and protect from light.

**Component:**

Primer pair (25 µM each) 180 µL

Probe (10 µM) 60 µL

**Detection of:**

Event- J3

**Probe & primer name:**

J3 F & J3 R

J3 P (FAM/TAMRA)

**Reference:**

1. "Testing for Foods Produced by Recombinant DNA Techniques" by Ministry of Health, Labour and Welfare, Japan.

**License:**

NIPPON GENE is permitted to manufacture and sell this product by the J.R. Simplot Company.



# GM バレイショ DNA J3 オリゴヌクレオチドセット

Code No. 318-08591 (1 Set)

**保存:**

Store at -20°C(遮光)

**内容:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 60 μL

**検出:**

Event- J3

**プローブ&プライマー名:**

J3 F & J3 R

J3 P (FAM/TAMRA)

**参考文献:**

1. 厚生労働省「安全性未審査の組換え DNA 技術応用食品の検査方法」

**License:**

ニッポンジーンは本品の製造および販売を J. R. Simplot より許諾されています。

本品は、試薬(試験研究用)として販売しているものです。  
医薬品の用途には使用しないでください。