



## GM Rice Detection (63Bt) Oligonucleotide Set 2

---

**Code No. 319-07541** (1 Set)

**Storage:**

Store at -20°C and protect from light.

**Component:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

**Detection of:**

63Bt

**Probe & primer name:**

T52-SF & OsNOS-R2

GM63-Taq (FAM/TAMRA)

**Usage example:**

Use 0.8 μL of this product (Primer pair) in each 25 μL reaction.

For detail on the PCR conditions, refer to the References.

**Reference:**

1. "Testing for Foods Produced by Recombinant DNA Techniques" by Ministry of Health, Labour and Welfare, Japan.

## GM コメ 63Bt コメ検出用オリゴヌクレオチドセット 2

**Code No. 319-07541** (1 Set)

**保存:**

Store at -20°C(遮光)

**内容:**

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

**検出:**

63Bt

**プローブ&プライマー名:**

T52-SF & OsNOS-R2

GM63-Taq (FAM/TAMRA)

**使用例:**

Bt コメ検出用のプライマー対として PCR 系 25 μL に本品 (Primer pair) を 0.8 μL 使用する。

PCR 条件の詳細については、参考文献を参照してください。

**参考文献:**

1. 厚生労働省「安全性未審査の組換え DNA 技術応用食品の検査方法」

本品は、試薬(試験研究用)として販売しているものです。  
医薬品の用途には使用しないでください。