



GM Rice Detection PLD (for endogenous gene) Oligonucleotide Set 2

Code No. 312-07531 (1 Set)

Storage:

Store at -20°C and protect from light.

Component:

Primer pair (25 µM each) 180 µL

Probe (10 µM) 180 µL

Detection of:

PLD (endogenous gene)

Probe & primer name:

PLD3959F & PLD4038R

PLD-P (FAM/TAMRA)

Usage example:

Use 0.8 µL of this product (Primer pair) in each 25 µL reaction.

For detail on the PCR conditions, refer to the References.

Reference:

1. "Testing for Foods Produced by Recombinant DNA Techniques" by Ministry of Health, Labour and Welfare, Japan.

コメ内在性 DNA PLD オリゴヌクレオチドセット 2

Code No. 312-07531 (1 Set)

保存:

Store at -20°C(遮光)

内容:

Primer pair (25 μM each) 180 μL

Probe (10 μM) 180 μL

検出:

PLD (コメ内在性遺伝子)

プローブ&プライマー名:

PLD3959F & PLD4038R

PLD-P (FAM/TAMRA)

使用例:

コメ陽性対照用プライマー対として PCR 系 25 μL に本品 (Primer pair) を 0.8 μL 使用する。

PCR 条件の詳細については、参考文献を参照してください。

参考文献:

1. 厚生労働省「安全性未審査の組換え DNA 技術応用食品の検査方法」

本品は、試薬(試験研究用)として販売しているものです。
医薬品の用途には使用しないでください。